




*i***STAND**[™]

infusion pump
stand



Height adjustable from
163cm to 307cm



Pole

Stainless steel resists staining, corrosion and rusting.

Base

Low center of gravity design and 12kg base offers improved stability.

58cm Diameter Base

Small base reduces stumbles and falls while patient ambulates.

P-1080-6 / iSTAND™

features

- Accommodates up to six infusion pumps.
- Bag Drop Prevention Device allows for controlled descent of pole.
- Stainless steel IV hook assembly features six hooks.
- Heavy duty 8cm sealed casters with precision ball bearings.
- Extremely durable powder coat finish base.
- Pole sleeve welded securely to the base for reinforced stability.



Accessories	Accessories
Patient Handle (Stainless) P-3301	Wire Basket (Stainless - Large)..... P-3201
Patient Handle (Black)..... P-3302-BLK For easy mobility.	Wire Basket (Black - Large) P-3202
Infusion Pump Holder (Multiple).....P-1080-C	Wire Basket (Stainless - Medium) P-3211
O2 Holder (E-size tanks) P-3601	Wire Basket (Stainless - Medium) P-3212

