



Double specialty back table decker



Maximize Space

Upper shelf creates additional space without compromising visibility.



Visibility

Easily arrange & organize instrument trays without the need for stacking.



Time

Improve setup efficiency & workflow during large surgical cases.



Drape

Create a sterile field using our unique **perfect fit™** drapes with transparent plastic backdrop.




double-decker + perfect fit drapes = *more usable OR space*



standard double-decker

- 427 122cm table
- 428 152cm table
- 429 183cm table

The OR Specific standard double-decker back table is ideal in supporting any surgical procedure that utilizes either one or two tiers, especially when space is limited.



adjustable double-decker

- 427-ADJ 122cm table
- 428-ADJ 152cm table
- 429-ADJ 183cm table

OR Specific's adjustable back tables accommodate variable staff heights ranging from 150cm to 198cm. The main surface of the table may be adjusted intraoperatively from a minimum height of 74cm to a maximum height of 100cm.

double-decker table features



T-bolts secure the upper tier once raised and positioned.

Table shown with upper shelf in storage position

Post hole for storing upper tier when not in use.

Two locking casters on front legs of table hold table stationary.

- **Adjustable Main Surface**
74cm min. height
100cm max. height
- **Adjustable Upper Shelf Positions**
30cm above main surface
38cm above main surface
- **Weight Capacity †**
68kg upper shelf
91kg main surface
91kg storage shelf
† weight uniformly distributed
- **Casters**
Four 9cm casters
Two front locking casters
- **Construction**
18 gauge stainless steel

perfect fit drapes

- Unique disposable drapes easily cover both the upper and main table surfaces.
- The clear plastic backdrop protects the sterile field while allowing light through to the main table surface.



international heavy duty drape

Three layer laminate system provides premium barrier protection, tear & puncture resistance, and a high capacity absorbent layer - 28 mil thick.

- 418-HDE-S 122 cm table 10 per case
- 419-HDE-S 152 cm table 10 per case
- 420-HDE-S 183 cm table 10 per case

Upper shelf creates additional space without compromising visibility or sterility.

Transparent plastic backdrop allows light to transfer to the main surface of the table for improved visibility.



- Three Velcro® tabs under the top tier remove any tenting effect and allow for complete access to the main surface.

